ENERGY DISPERSIVE X-RAY FLORESCENCE (EDXRF)

	<u>ORESCENCE (EDXRF) BENCHTOP SPECTROMETER OR</u> QUIVALENT:	
<u>S</u>	pecifications:	
•	Element range Carbon (C) to Americium (Am) 15 Watt version with high current (3mA): used for higher sample throughput with improved and extended light element capabilities Carbon (C) to Americium (Am) Can handle 10 samples at a time Can handle sample type few milligram to larger bulk samples 10-position removable sample changer Spectrometer can accommodate samples up to 52 mm (2 inch) in diameter Spinner is included for better accuracy of air filter analysis Large sample mode for analyzing bigger samples up to 10 cm in height	
	Epsilon Software PTIONAL ACCESSORIES:	
	High-resolution silicon drift detector (SDD30), typically 135 eV @ Mn-Ka	01
	Elemental Screening OMNIAN software WROXI calibration standard	01 01
<u>P</u>]	ICTURE OF EPSILON 4XRF	•

- **d.** Certificate from the manufacturer that all the items are brand new & of good quality.
- **e.** On-site inspection & training.
- **f.** Quotation must be accompanied with full technical literature.