

GC-MS SPARES AND CONSUMABLE

Sr. No.	Description/ Specifications/	Part No	Qty. Required
1	6900MS Kit	GCMSPMKIT	05 kit
2	Rotary Vane Pump Kit of Edwards	GCRPPMKIT	05kit
3	Calibration Sample	GCMSCAL	05 kit
4	μ-ECD Assembly for YL65Q0GC Kit	YGEDPMKIT	02kit
5	Set of spacers for 40 position tray to allow use of 10ml vials	HT200H-2010	01 set
6	20ml Headspace vials clear, round bottom, 100ea/pkt	1.29.942	20pk
7	20ml vials, Crimp Cap with Sil/PTFE septa, 100/pack.	1.29.940	20pkt
8	10ml vials, Crimp Cap with Sil/PTFE septa, 100/pack.	1.29.940	20pkt
9	Plunger only for 2.5ml Syringe	2005622690	05ea
10	Glass Barrel + integral needle for 2.5ml syringe	2C05622700	05ea
11	Hand Crimper, 5-100ml (20mm seal)	33280-U	02ea
12	Decapper for 20mm seals, EA	33183	02ea
13	Auto sampler syringe. 10 pi	1602633490	10ea
14	Cans for wash/waste vials. Pack/20	1.29965	40pkt
15	Replacement Septa for Wash/Waste vials 20ea/pk	1.29.936	40pkt
16	Big Universal Trap, 1/8" ft, Helium	RMSH-2	01ea
17	10 μl syringe, 6/pk, manual type	1602(533480	05pkt
18	Bleed/Temp Optimized Inlt Septa 11mm. 50pk (upto 400 °C)	518.- -475 7	05pkt
19	Liner, Ultra Inert, split less, single taper, glass wool	5190-2293	05ea
20	Liner. Ultra Inert, universal, low pressure drop, glass wool	5190-2295	10ea
21	O-Ring. inlet liner. non-stick. 10/pk	1321630034	10pkt
22	Ferrule for MSD connection, 10ea/pk (0.32mm ID capillary column)	SGE-072654	02pkt
23	Ferrule for MSD connection, 10ea/pk (0.25mm ID capillary column)	SGE-072663	05pkt
24	NM-Hi Plus, 250(Hydrogen generator, flow 250ml/min), 2.20-240V/50-60Hz	2103923110	01ea
25	Capillary Column for CMME analysis (with application note)	CAPCOLUM	01ea
26	DB-lms Ultra Inert 30m, 0.25mm, 0.25pm (with application note)	122-0132UI	01ea

GC-MS SPARES AND CONSUMABLE

27	DB-5M5 30m, 0.25mm, 0.25µm (with application note)	122-5532	01ea
28	DB-624 UI 30m, 0.25mm, 1.40µm (with application note)	122-1334UI	01ea
29	DB-17MS 30m, 0.25mm, 0.25µm (with application note)	122-4732	01ea
30	DB-WAX Ultra Inert 30m, 0.25mm, 0.25µm (with application note)	122-7032UI	01ea
31	Pulse Discharge Detector for YL6500 GC, 220V, HID mode for capillary Direct Control and acquisition	1351620060	01ea
32	NPDOTID) for 6500GC	6501011071	01ea
	Make: YL Instrument Co. Korea		