Sr. No.	Item	Quantity
1.	Dosimat With Titration Stand	03No.
	Complete with Exchange units of 20 and 50 ml burette cylinder volumes	
	Capable of dispensing in different modes	**
	Capable of controlling the volume	
	• 1 liter reagent bottle of amber glass with two burette tips	
	<ul> <li>Anti-diffusion tips to prevent outside liquid from diffusing inward</li> </ul>	
	Display LCD in variable contrast	
	Magnetic stirrer with titration stand and controllable speed	
	• Power supply: 220V, 50Hz.	
	Accessories	
	Connecting cables and tubes etc	0.7
	Spares	each
	Internal fixed electronic cards, Cables, Connecting tubes, Anti-	each
	diffusion tips, dosing tips open	,
	Make: Metrohm, Thermo or Mettler toledo.	
	> (Technical Literature must be provided with quotation)	