

Sr. No.	Item	Quantity
1.	<p><b><u>Dosimat With Titration Stand</u></b></p> <ul style="list-style-type: none"> <li>• Complete with Exchange units of 20 and 50 ml burette cylinder volumes</li> <li>• Capable of dispensing in different modes</li> <li>• Capable of controlling the volume</li> <li>• 1 liter reagent bottle of amber glass with two burette tips</li> <li>• Anti-diffusion tips to prevent outside liquid from diffusing inward</li> <li>• Display LCD in variable contrast</li> <li>• Magnetic stirrer with titration stand and controllable speed</li> <li>• Power supply: 220V, 50Hz.</li> </ul> <p><b><u>Accessories</u></b> Connecting cables and tubes etc</p> <p><b><u>Spares</u></b> Internal fixed electronic cards, Cables, Connecting tubes, Anti-diffusion tips, dosing tips open</p> <p><b>Make:</b> Metrohm, Thermo or Mettler toledo.</p> <p>➤ (Technical Literature must be provided with quotation)</p>	<p>03No.</p> <p>02 each</p>